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
H&E Electronic Journal

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^{*}Registered trademark; "Vitravene" by ISIS Pharmaceuticals, "Macugen" by Pfizer, "Kynamro" by Sanofi Genzyme, "Exondys 51" by Sarepta Therapeutics, "Spinraza" by Biogen, "Defitelio" by Jazz Pharmaceuticals, "Tegsedi" by Ionis Pharmaceuticals, "Onpatro" by Alnylam Pharmaceuticals, "Hepilisav-B" by Dynavax Technologies, "Waylivra" by Ionis Pharmaceuticals, "GIVLAAR" by Alnylam Pharmaceuticals, "VYONDYS 53" by Sarepta Therapeutics, "Viltepso" by Nippon Shinyaku Co. Ltd.

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
Fig.1 Market size and forecast of global nucleic acid drug between 2018 and 2024

Source: Evaluate Pharma/ Boston Consulting Group/ China Galaxy Securities

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
Ref.; Pharmaceutical Technology, February 2020, p 30-34, <https://www.pharmtech.com/view/designing-commercial-scale-oligonucleotide-synthesis>

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
Ref.; Pharmaceutical Technology, February 2020, p 30-34
<https://www.pharmtech.com/view/designing-commercial-scale-oligonucleotide-synthesis>

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The elution profile for the TSKgel SuperQ-5PW (20) resin went off scale at 254 nm. The enlarged image of the main oligonucleotide peak, when overlaid with a histogram showing the AIEC HPLC results for fraction purity, highlights the chromatographic separation of the full length oligonucleotide that could not be determined by examining the elution profile alone. Note: This figure highlights just the area in full length oligonucleotide elutes.

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Resin	Crude Loaded	Crude Purity (AEX)	Crude Purity (RPC)	Total Recovered (OD 260)	% Full Length Oligo (AEX)	% Full Length Oligo (RPC)	% Yield of FLO (AEX)	% Yield of FLO (RPC)
TSKgel SuperQ-5PW (20)	37.1 mg/mL	66.5%	62.96%	118.8 mg	97.3%	86.3%	74.3%	69.6%
TOYOPEARL SuperQ-650S	45.0 mg/ mL	66.5%	62.96%	142.7 mg	95.6%	82.8%	74.1%	67.7%
TOYOPEARL GigaCap Q-650S	30.0 mg/ mL	66.5%	62.96%	95.3 mg	96.5%	85.5%	75.0%	70.0%

Sample purity was determined by pooling the highest purity fractions (>85% full length oligonucleotide - FLO) from each column and analyzing them on a TSKgel DNA-STAT HPLC column. The pooled fractions were also analyzed using a Waters Acquity OST C18 1.7 μ m (2.1 mm ID x 10 cm) RP UPLC column. Recovery was determined by comparing the amount of full length oligonucleotide present in the crude sample loaded onto the column with the amount of full length oligonucleotide present in the pool of high purity fractions.


- TSKgel SuperQ-5PW (20) showed slightly better results on yield and recovery than TOYOPEARL GigaCap Q-650S.
- TOYOPEARL SuperQ-650S showed higher dynamic binding capacity, while the yield and recovery were slightly low among the three anion-exchangers examined.
- TOYOPEARL GigaCap Q-650S showed 30 % lower volume for elution pool of the sample compared with other anion-exchangers..

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

Ref.; Tosoh Bioscience LLC Application Note No. 117

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Hybridon Inc., Patent WO 98/01288 (1996)
Biogen MA Inc., Patent WO 2017/218454 A1 (2017)