TSK-GEL® Super-Phenyl Products

Part Numbers: 18277, 4.6mm ID X 5cm, 2μm

 $18278, 4.6 mm \ ID \ X \ 10 cm, 2 \mu m$ 20017, 2.0mm ID X 5cm, 2μm 20018, 2.0mm ID X 10cm, 2μm 18207. Guard filter (3/pk) 18206, Guard holder

This sheet contains the recommended operating conditions and the specifications for TSK-GEL Super-Phenyl columns and guard filters Installation instructions and column care information are described in a separate Instruction Manual.

OPERATING CONDITIONS

1. Shipping Solvent: 50% Acetonitrile - 50% Water

2. Max. Flow Rate: 4.0 mL/min (4.6mm ID) 0.25 mL/min (2.0mm ID)

When a buffer with high viscosity is used, the maximum flow rate may have to be reduced so it doesn't exceed the

maximum pressure drop. When changing solvents, use a flow rate equal to 25% of the maximum flow.

Standard Flow Rate: 1.0 - 2.5 mL/min (4.6mm ID)

0.15 - 0.20 mL/min (2.0mm ID) 30.0 MPa (4.6mm ID) Max. Pressure:

8.0 MPa (2.0mm ID x 5cm)

15.0 MPa (2.0mm ID x 10cm)

pH Range: 2.0 - 7.5

Organic Conc. Range: No limitation

Temperature: 10 - 50°C. Reduce flow rate when operating below 10°C.

Cleaning Solvents: (1) High conc. solvent containing organic modifiers*

(2) Mixture of organic acids and high conc. organic modifiers*

*Acetonitrile and methanol are recommended as a modifier.

Store the column in the shipping solvent if it will not be used within three days. Prevent air from entering the 9. Storage:

column. For overnight storage flush the column with mobile phase at 0.2mL/min.

An on-line filter (0.2-0.5 $\mu m)$ between pump and injection valve is recommended. Guard filters prevent the column Column Protection:

from a contamination of strongly adsorbed solutes. As a general rule, guard filters should be replaced after 30-40

sample injections or when peaks become excessively wide.

SPECIFICATIONS

The performance of TSK-GEL Super-Phenyl column is tested under the conditions described in the data sheet. All columns have passed the following quality control specifications:

≥ 8,000 (4.6mm ID x 5cm) 1. Number of Theoretical Plates (N):

≥ 16,000 (4.6mm ID x 10cm) \geq 3,000 (2.0mm ID x 5cm)

≥ 6,000 (2.0mm ID x 10cm)

2. Asymmetry Factor (AF): 0.8 - 1.6 (4.6mm ID)

0.8 - 1.8 (2.0mm ID x 5cm) 0.7 – 1.6 (2.0mm ID x 10cm)