TSK-GEL® ODS-140HTP Products

21927, TSKgel ODS-140HTP, 2.1mm ID x 5cm, 2.3 μ m 21928, TSKgel ODS-140HTP, 2.1mm ID x 10cm, Part Numbers:

 $2.3 \mu m$

This sheet contains the recommended operating conditions and the specifications for TSK-GEL ODS-140HTP columns. Column care information is described in a separate Instruction Manual.

OPERATING CONDITIONS

3. pH Range:

1. Shipping Solvent: Acetonitrile 2. Max. Pressure: 60MPa

4. Organic Conc. Range: 2 - 100%

5. Temperature: 10 - 50°C. Reduce flow rate when operating below 10°C.

6. Cleaning Solvents: (1) High conc. solvent containing organic modifiers*

2.0 - 7.5

(2) Mixture of organic acids and high conc. organic modifiers*

*Acetonitrile and methanol are recommended as a modifier.

7. Storage: Store the column in the shipping solvent if it will not be used within three days. Prevent air from

entering the column. For overnight storage flush at 0.1mL/min.

SPECIFICATIONS

The performance of TSK-GEL ODS-140HTP columns are tested under the conditions described in the data sheet. All columns have passed the following quality control specifications:

1. Number of Theoretical ≥ 7000 (2.1mm ID x 5cm) Plates (N): > 14,000 (2.1mm ID x 10cm)

2. Asymmetry Factor (AF): 0.9 - 1.3 (2.1mm ID)